

Measuring
people's
safety

second
breath 

FASTER MORE FLEXIBLE EASIER


**TEST EQUIPMENT
FOR THE RESPIRATORY
PROTECTIVE DEVICES**



About company

SECOND BREATH OÜ is a research and manufacturing company based in Estonia whose subject of activity is innovative lab equipment and software for testing and quality control of all types of RPE, monitoring the physiological parameters of human breathing.

Our products are recommended for:



About company

**5+
YEARS**

of experience in own
production of test equipment
for RPE

**14 SERIAL
TEST
EQUIPMENT**

in the product line

**25+
YEARS**

of experience of the
company's employees
in RPE research
and development

**20+
EMPLOYEES**

in our team including engineers,
designers, programmers,
brand manager, customer service
and sales managers

**70%
OF EXPERTS**

have a scientific
degree
in engineering
and IT

6 NEW PRODUCTS

planned for release in 2024

How we manufacture test equipment



Market,
customer
requests
analysis



R&D
work



Components
selection from
the world's best
manufacturers



Own unique
hardware and
software parts
development



Testing and
internal quality
control



Certification
and obtaining
a patent
for innovative
technologies



Market
launch and
direct sales



SECOND BREATH OÜ and mining industry

20+ YEARS

our main engineer has experience in research and development SCSR



Metabolic breathing simulator OXY ROBOT

- for testing SCSR
- for incoming and periodic inspection
- for annual monitoring SCSR program



is the first equipment we produced

Mobile breath recorder

for testing SCSR with volunteer person



is unique device which has no analog in the world

Recommended for:

- developers and manufacturers SCSRs
- mining companies and mine rescue teams
- scientific and technology research and control organization and lab

Test equipment for breathing apparatuses and life support systems

OXY ROBOT

metabolic
breathing
simulator



OXY ROBOT

pressure breathing simulator



MOBILE BREATH RECORDER



PULSAR

breathing machine



OLYMPUS

pressure breathing
machine



DUMMY HUMAN HEAD




Test equipment for filtering devices

BOREAS
measuring
breathing resistance



PULSAR
breathing machine



OLYMPUS
pressure breathing
machine



DIOXIDE
measuring CO₂
content in the
«dead space»



MONOXIDE
measuring CO
content in inhaled air



DOLOMITE
measuring filters resistance
to dolomite dust clogging



SILICA DUST
measuring filters resistance
to silica dust clogging



DYNAMICS
evaluation of
gas, vapor and
combined
filters
capacity



SALT MIST
measuring total inward leakage (TIL)
with using NaCl aerosol



ELEGAS
measuring total inward leakage (TIL)
with using SF₆



DUMMY HUMAN HEAD
testing different facepieces



SECOND BREATH OÜ test equipment features:

WE

- test equipment comes to working mode fast
- test equipment have a high level of process automation and do not require constant staff involvement

WE

- test parameters can be flexibly changed* in the test equipment software via an application on a personal computer
- * within technical possibilities*

WE

- test equipment have a simple and convenient interface to control

WE

- test equipment have lightened design

It's FASTER

It's MORE FLEXIBLE

It's EASIER

It has ERGONOMIC DESIGN

YOU

- can carry out more tests per day
- increase your productivity
- do other important work task at the same time

YOU

- not have to change all or part of the test equipment if industrial standards are updated or new one release, or if you carry out research tests


YOU

- get minimum manual switching
- can quickly and easily master the control of the test equipment and involve even a laboratory intern in the work

YOU

- get attractive laboratory and additional inspiration to design RPE

About service




**Below you will take a closer look at our equipment,
but first, remember a few icons that will help you to navigate**



for testing
breathing
apparatuses




for testing
filtering devices



work in a
climate chamber



work in
conditions
of low or high
atmospheric
pressure



testing with
volunteer
person



Human breathing simulation

Flagship product

OXY ROBOT metabolic breathing simulator



recreate human inhale and exhale with necessary humidity and temperature, the exchange of gasses: O_2 consumption, CO_2 release; reproduce breathing according to prerecorded breathing curve of a real person at different loads.
Can be completed with a pumped or heated dummy human head or tee.

APPLICATION:

SCBA, Airline BA, SCSR on chemically bonded or compressed oxygen, EBA, life support systems, CCR, SCUBA

RELEVANT STANDARDS*:

ISO 16900-5:2016(E), EN 401, EN 14143 and others*




Watch clip


OXY ROBOT pressure breathing simulator




suitable for tests in conditions of ambient pressure difference (6 atm).
Design with cylinders



Controlled
via touch-screen
display



Testing breathing apparatuses with persons



MOBILE BREATH RECORDER test equipment




Complements the testing program conducted with Metabolic breathing simulator OXY ROBOT

a portable wireless device and software that is connected to a real RPE and records the parameters of a person's breathing in it under various loads and climatic conditions (up -40 to 40 C°).

the load during testing process. You can create an exercise program for testing and find out the breathing parameters at the time of these loads from the report.

smartphone and application for choosing the loads and getting feedback about own breathing

MEASUREMENT PARAMETERS:



Watch clip

tested breathing apparatus

portable measuring unit — unique system connects with RPE and record breathing parameters

volunteer person



If you interest in this products, send your request on info@second-breath.ee to get the datasheet

Lung ventilation simulation




Artificial lungs

PULSAR breathing machine



create a pulsating air flow that imitates human breathing with the required value of lung ventilation. Can be completed with a pumped or heated dummy human head.



Controlled by personal computer or smartphone via Bluetooth-channel

suitable for tests at normal atmospheric pressure. Design with bellows pumps

OLYMPUS pressure breathing machine




suitable for tests in conditions of ambient pressure difference. Design with cylinders

Measuring breathing resistance

Test device BOREAS



Controlled
via touch-screen display



creates a constant air flow through the RPE in the directions of «exhalation» and «inhalation» and measures the air pressure difference between the surrounding atmosphere and the mask space of the RPE, calculates the resistance to air flow.

APPLICATION:

RPE with full-face and half-mask; anti-aerosol, anti-gas and combined RPE with an insulating front part; filtering self-rescuer; RPE from radioactive substances.

RELEVANT STANDARDS*:

EN 136, EN 140, EN 143, EN 149, EN 405, EN 1827 and others...



Watch clip

Evaluation of filters capacity using test agents



Test equipment DYNAMICS



Watch clip

prepares the gas-vapor-air mixture in accordance with the conditions for using the tested filter - the required temperature, relative humidity, concentration of the test agents and passes this mixture through the filter until the test agent is detected in the air sample at the outlet of the filter (breakthrough concentration) or until the filter protection time, which is set by the RPE manufacturer, expires.

APPLICATION:
gas, vapor and
combined filters



RELEVANT STANDARDS*:
EN 136, EN 140, EN 143, EN 149, EN 405, EN 1827
and others...

Test different types of filters on the same test equipment. Configuration options for one test equipment design:

1) Dynamics G: up 2 to 8* test gases:





	Cl ₂	chlorine		C ₂ H ₆ O	dimethyl ether
	NO ₂	nitrogen dioxide		C ₄ H ₁₀	isobutane
	SO ₂	sulfur dioxide		H ₂ S	hydrogen sulfide
	NH ₃	ammonia		NO	nitric oxide

2) Dynamics G2: 2* gases:

	HCl	hydrogen chloride
	AsH ₃	arsenic hydrogen

* the number of test agents
is determined by the customer

3) Dynamics V: up 1 to 4* vapors

	C ₆ H ₁₂	cyclohexane
	HCN	hydrogen cyanide
	C ₃ H ₄ O	acrolein
	C ₆ H ₆	benzene

4) Dynamics Hg

	Hg	hydrargyrum
---	----	-------------



Controlled via personal computer (laptop)

Measuring CO₂ and CO content in inhaled air



developed on OXY ROBOT breathing simulator hardware and software platform

DIOXIDE test equipment




creates a pulsating air flow that imitates human breathing with the necessary parameters of breathing frequency and depth and volume concentration of CO₂ on «exhalation» and measures CO₂ content in a mask “dead space”.

APPLICATION:

full-face masks, half-masks, quarter-face masks, masks with helmet or hood

RELEVANT STANDARDS*:

EN 1827, EN 12941, EN 136, EN 149, EN 403, EN 404, EN 405 and others...



Controlled via touch-screen display



Watch clip



Watch clip

creates a pulsating air flow that imitates human breathing with necessary parameters of breathing frequency and depth, temperature and relative humidity, creates a test atmosphere with necessary conditions — temperature, moisture content, volume of carbon monoxide (CO) and measures CO content that has passed through the filter of tested RPE in the required time.

APPLICATION:

filtering self-rescuers

RELEVANT STANDARDS*:

EN 403, EN 404 and others...

MONOXIDE test equipment



Measuring dust clogging resistance

developed on OXY ROBOT breathing simulator hardware and software platform



DOLOMITE test equipment acc. to EN standards



Using test agent is dolomite dust

test equipment creates a special dust atmosphere and constant air flow in the test chamber for the tested respirator and the breathing machine simulates human breathing with necessary parameters of lung ventilation, humidity and temperature of exhalation. As a result of the test respirator resistance to constant airflow after clogging is determined. This time is compared with requirements of industry standards.



Controlled via personal computer (laptop)



Watch clip



Watch clip

SILICA DUST test equipment acc. to NIOSH standard



Using test agent is silica dust

APPLICATION:

respirators with anti aerosol
and combined filters

RELEVANT STANDARDS*:

EN 149, EN 143, EN 13274-8, EN 405,
EN 1827, EN 12941, EN 12942, EN 14387
and others...

APPLICATION:

HEPA filters designed for powered air-purifying respirators
(PAPRs)

RELEVANT STANDARDS:


procedure № RCT-APR-STR-0025 NIOSH

Measuring total inward leakage (TIL)

test RPE with persons



ELEGAS test equipment



Using test agent
is sulfur hexafluoride (SF6)



Watch clip

test equipment creates a special atmosphere with a test agent and determines the leakage of it under a mask bypassing a filter when a volunteer person is doing exercises in the tested RPE on the treadmill into the test chamber.

APPLICATION:

full face & half masks with aerosol filter, half masks with valves and non-removable gas or combined filter, filtering self-rescuers

RELEVANT STANDARDS*:

EN 149, EN 13274-8, EN 136,
EN 1827, EN 12941, EN 12942
and others...



Controlled via personal computer

SALT MIST test equipment




Using test agent
is sodium chloride aerosol (NaCl)



Watch clip

Dummy human head

sold only as part of Oxy Robot, Pulsar, Olympus testing equipment




**Small
ISO size**



**Medium ISO
size**



**Large
ISO size**




+ customization size is available
too according to your 3D model

APPLICATION:

1. for testing various facepieces - full-face, half-face, with a valve, in the form of a hood
2. can be made in a heated version for climatic chambers (without inflation)

How quickly SECOND BREATH OÜ test equipment reach* operating mode



* Time to reach the operating mode is the time interval after power supply is applied to the equipment, after which the specified technical characteristics required for testing must be achieved.

Measuring
people's
safety



- 📍 Address: Valge 13, 11415, Tallinn, Estonia
- ☎ Phone & business WA: +37 269 807 09
- ✉ E-mail: info@second-breath.ee
- 🌐 Web: www.second-breath.ee

STAY IN THE LOOP!

- 🌐 LinkedIn: Second Breath OÜ
- 📺 Youtube: [@secondbreath5264](https://www.youtube.com/@secondbreath5264)

VIDEO PRESENTATION

